

Wireless LED Curing Light LED.B

Woodpecker · Model LED.B — Curing Lights

Lightweight 145 g wireless curing light with a replaceable lithium battery. Light intensity stays constant as the battery drains, so curing depth is unaffected by remaining charge, and standby time reaches roughly forty days.

Key features

- Wireless with replaceable lithium battery
- Timer settings of 5 s, 10 s, 15 s and 20 s
- 400+ cures per full charge in 10 s mode
- Low standby consumption, ~40 day standby
- Constant light intensity throughout discharge

Technical specifications

Light output	850–1000 mW/cm ²
Wavelength	420–480 nm (blue light)
Power input	AC 100–240 V, 50/60 Hz
Consumption	≤8 W
Battery	3.7 V / 2000 mAh rechargeable lithium
Dimensions	31 × 34 × 260 mm
Net weight	145 g
Protection class	Class II, IPX0

Applications

- Composite restorations
- Adhesive bonding
- Fissure sealants

Supplied, installed and serviced across Nepal by Serving Nation Enterprises. Head Office: Putalisadak, Kathmandu 44600 · Branch Office: Chitwan, Nepal. Specifications subject to change without notice.