

Integrated Dental Chair Unit TR-T5

Tree Dental · Model TR-T5 — Dental Chairs

Air-controlled integrated dental unit with 24V noiseless DC motor, thermostatic water supply, foldable biarticular headrest and a choice of two LED sensor operating lamps. Side box and patient chair share one integrated design for a compact footprint.

Key features

- 3 memory positions with air-brake instrument tray
- 24V DC noiseless motor and air control system
- Automatic thermostatic water supply and rinse control
- Double rotatable luxury armrests for easy patient access
- Built-in purified water system with 500 ml bottle
- Built-in or external floor box options

Technical specifications

Operating lamp	LED sensor lamp, 6500–26000 lux, 4300 K
Alternate lamp	6-bulb LED, 8000–30000 lux, 5000 K, 9 W
Handpiece tubing	3 pcs, 4-hole or 2-hole
Syringes	2 × 3-way syringe (hot/cold)
Spittoon	Removable glass, rotates 90°
Included	LED film viewer, dentist stool, multifunction foot controller

Applications

- Startup dental clinics
- Community dental centres
- Dental camps

Supplied, installed and serviced across Nepal by Serving Nation Enterprises. Head Office: Putalisadak, Kathmandu 44600 · Branch Office: Chitwan, Nepal. Specifications subject to change without notice.