

European Dental Chair Unit TR-E1

Tree Dental · Model TR-E1 — Dental Chairs

Motor-driven aluminium base dental chair with two memory positions (spit and recovery), wide clear operating board, 45°-rotatable European unit box and a double water bottle system (600 ml / 1000 ml) that also accepts disinfecting solution for line cleaning.

Key features

- Motor driving aluminium base with memory positions
- Swing-mount delivery system, down- or top-mounted tray
- Double-joint headrest and barrier-free rotating armrests
- Pressure-stabilising valve keeps handpiece speed constant
- Automatic water heating, cup filler and spittoon flush
- Assistant rack with strong/weak suction, rotates 120°

Technical specifications

Delivery system	Swing mount, electric motor driving
Handpiece tubing	2 × high speed, 1 × low speed
Syringes	2 × 3-way syringe
Water bottles	600 ml + 1000 ml clean water system
Spittoon	Deluxe ceramic turntable, rotates 90° inward
Included	Operating lamp, 24V film viewer, ergonomic dentist stool

Applications

- General dentistry
- Restorative and endodontic treatment
- Dental teaching clinics

Supplied, installed and serviced across Nepal by Serving Nation Enterprises. Head Office: Putalisadak, Kathmandu 44600 · Branch Office: Chitwan, Nepal. Specifications subject to change without notice.