

North American Style Dental Chair

TR-S105

Tree Dental · Model TR-S105 — Dental Chairs

Air-controlled unit with imported DC motor, 180° swing-mount delivery system, movable side box rotating 60° and a three-stage water filtration system. Microfiber leather upholstery with lumbar pad support for long procedures.

Key features

- 180° swing mount delivery system
- Movable side box rotating at least 60°
- Assistant element rotates 90°
- Three-stage water filtration system
- Six-bulb induction oral lamp
- Up to 9 memory positions

Technical specifications

Upholstery	Microfiber leather with lumbar pad
Cuspidor	Turnable ceramic
Handpiece tubing	3 pcs, 4-hole or 2-hole
Syringes	2 × 3-way syringe (cold and hot)
Water heater	Low-voltage AC
Included	Film viewer, memory foam dentist stool, floor box

Applications

- Dental clinics
- Orthodontic practices
- Hospital dental departments

Supplied, installed and serviced across Nepal by Serving Nation Enterprises. Head Office: Putalisadak, Kathmandu 44600 · Branch Office: Chitwan, Nepal. Specifications subject to change without notice.