

Noiseless Oil-Free Air Compressor

TR-OAC01

Tree Dental · Model TR-OAC01 — Dental Compressors

Complete compressor family from a 600 W single-chair unit to a 6000 W 25-chair system. The optional condensing air dryer freeze-dries compressor air, then a molecular-sieve dryer and air filter purify it, with an intelligent automatic drainage system that removes the need for manual draining.

Key features

- Oil-free operation — no oil contamination in the air line
- Low noise level from 53 dB
- Optional condensing air dryer with molecular sieve
- Intelligent automatic drainage system
- Air filter keeps supply pure and motor quiet
- Models for 1, 2, 3, 4, 5, 6, 7, 8, 10, 15, 20 and 25 chairs

Technical specifications

TR-OAC01-1	1 chair · 600 W · 118 L/min · 24 L tank · 53 dB
TR-OAC01-2	2 chairs · 800 W · 155 L/min · 40 L tank · 53 dB
TR-OAC01-3	3 chairs · 1200 W · 236 L/min · 60 L tank · 53 dB
TR-OAC01-4	4 chairs · 1600 W · 310 L/min · 60 L tank · 54 dB
TR-OAC01-5	5 chairs · 1800 W · 354 L/min · 90 L tank · 55 dB
TR-OAC01-6	6 chairs · 2400 W · 465 L/min · 90 L tank · 55 dB
TR-OAC01-8	8 chairs · 3200 W · 620 L/min · 120 L tank · 56 dB
TR-OAC01-10	10 chairs · 4000 W · 775 L/min · 150 L tank · 58 dB
TR-OAC01-20	20 chairs · 4800 W · 960 L/min · 200 L tank · 380 V
TR-OAC01-25	25 chairs · 6000 W · 1200 L/min · 240 L tank · 380 V
Voltage	110 V / 220 V, 50/60 Hz (380 V on large models)

Applications

- Dental chair air supply
- Dental laboratories
- Hospital dental departments

Supplied, installed and serviced across Nepal by Serving Nation Enterprises. Head Office: Putalisadak, Kathmandu 44600 · Branch Office: Chitwan, Nepal. Specifications subject to change without notice.