

UV Disinfection Cabinet TR-256N

Tree Dental · Model TR-256N — Sterilizers & Disinfection

Tempered-glass UV storage unit that disinfects and stores instruments between procedures. The digital version adds timing control and an ozone disinfection function; both versions can be wall-mounted to save floor space.

Key features

- 253.7 nm ultraviolet disinfection
- Digital version with timer and ozone function
- Tempered glass see-through doors
- Wall-mounting holes for space saving
- Micro switch cuts UV when the door opens
- 10 metal pallets per door (25 × 15 cm)

Technical specifications

Capacity	27 L single door · 27 L × 2 double door
Power	15 W
Voltage	220 V, 50 Hz
UV wavelength	253.7 nm
Single door size	310 × 230 × 480 mm
Double door size	640 × 230 × 480 mm
Weight	≈10 kg single · ≈20 kg double

Applications

- Dental instrument disinfection
- Sterile storage
- Clinic and lab infection control

Supplied, installed and serviced across Nepal by Serving Nation Enterprises. Head Office: Putalisadak, Kathmandu 44600 · Branch Office: Chitwan, Nepal. Specifications subject to change without notice.