

Dental Suction Machine TR-YP606A1

Tree Dental · Model TR-YP606A1 — Suction & Vacuum Systems

Compact central suction unit available in 550 W (2 chairs) and 750 W (3 chairs) versions, running in intermittent or continuous operation mode for clinics that need reliable shared aspiration.

Key features

- Two-chair and three-chair versions
- Intermittent and continuous operation modes
- 32 mm inlet and drain pipe connections
- Stable vacuum for simultaneous chairs
- Compact 38 × 34 × 60 cm footprint

Technical specifications

2-chair model	550 W · ≥ 1.1 kPa · 400 L/min · 42 kg
3-chair model	750 W · ≥ 1.2 kPa · 600 L/min · 43 kg
Voltage	220 V 50 Hz / 110 V 60 Hz
Pipe diameter	32 mm inlet and drain
Product size	38 × 34 × 60 cm
Ambient conditions	+5 °C to +40 °C, RH $\leq$ 80%

Applications

- Shared clinic suction
- Oral surgery rooms
- Hospital dental units

Supplied, installed and serviced across Nepal by Serving Nation Enterprises. Head Office: Putalisadak, Kathmandu 44600 · Branch Office: Chitwan, Nepal. Specifications subject to change without notice.